# Powertube®

DNV

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**Seite 3**

**Mode d'emploi**

**page 18**

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**pagina 33**

**Instruction manual**

**page 48**

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---

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# Instruction manual

---

***This instruction manual is part of the product. It must be kept handy near the device at all times. Strictly following this manual is essential to ensure the proper use and correct handling of the device, and for the safety of the patient and the user.***


---

## Table of contents Page


1.	What does the device look like?	4
2.	Control panel with functions	5
3.	Description of the device	6
4.	Quick guide to starting up	6
5.	Battery and battery holder	7
6.	Technical specifications	8
7.	How to use the device	9
8.	Indications	10
9.	Contraindications	10
10.	Safety instructions for use	10
11.	Care, maintenance, and storage	11
12.	Errors and faults, function test	12
13.	Accessories	12
14.	Guarantee	12
15.	Safety check	13
16.	EC Declaration of Conformity	14

# 1. What does the device look like?

## 2. Control panel with functions

1.  Switch on the device (press briefly)
2.  Select minutes (1 - 7 minutes) Display
3.  Start therapy (press briefly) (STEP manual selection)
4.  Switch off the device (long press or automatically after 30 seconds)

## 3. Description of the device

### Purpose

POWERTUBE is used for electronic nerve stimulation. It is used for treating pain at home.

This device is not a replacement for treatment or guidance from a doctor, but it may alleviate pain or provide relief from pain. Instructions for use are given in detail, in the section “How to use the device”.

### Functional description

The TENS (transcutaneous electrical nerve stimulation) device POWERTUBE is an easy-to-use, battery-operated, electronic device for nerve stimulation.

Electrical impulses are used to directly stimulate the corresponding nerve points. The device can be used by persons who are not medically trained.

This self-treatment is safe and has no side effects. The attached instruction manual provides an easy-to-understand description of the device and makes self-treatment simple. The device is designed for use by persons who have no medical training. The instruction manual shows the possible uses and their duration.

## 4. Quick guide to starting up

### Before first use

The device should be checked for transport damage as soon as it is received and subjected to a functional test. If damage has occurred during transport, a complaint must be filed with the carrier immediately and a damage report must be made to secure the claim for compensation.

In addition to the name and address of the recipient, this must also include the date of receipt, type and serial number of the delivered device, and description of the damage. The original packaging of the device should be kept during the warranty period so that the device can be returned in its original packaging if necessary.

Before delivery, each device is tested by the manufacturer for its function and safety. To ensure that the device continues to function safely after transport and installation by the user, these instructions should be followed:

---

***The user should not put the device into operation until a functional test has been carried out (see section 12).***

---

## 5. Battery and battery holder

Use alkaline batteries and observe the polarity (plus/minus). These batteries ensure that therapy can be completed even if the “low battery” indicator lights up when the battery is almost empty. The POWERTUBE devices are designed for this type of battery. **DO NOT** use lithium batteries.

**Do NOT** use rechargeable batteries because the device will not function properly with these batteries.

## 6. Technical specifications

Display	Therapy time in minutes, 7 LEDs, Therapy steps 1 - 3, replace LED battery.
Error message	Interval sound: Hand electrode or contact surfaces not applied correctly.  Continuous sound: Short circuit between hand electrodes and the POWERTUBE body.
Battery	9 Volt Alkaline Block 6LR61
Output voltage	AC max. 9 Volt pp and min. 6 Volt pp
Therapy time	3 - 21 minutes
Programme sequence	3 steps
Dimensions	Length 18 cm, diameter 3.1 cm
Weight in soft bag	≈ 750gr
Device classification	Type BF 
Classification	Ila Device
IP protection classes	Device IP40 Protection against penetration of solid foreign bodies with a $\varnothing > 1$ mm, not protected against water.  Soft bag: IPX2 Protection against water falling at an angle (dripping water), 15° compared to normal operating position.
Operation	<ul style="list-style-type: none"> <li>• 5°C to 40°C</li> <li>• Humidity 15% - 93% without condensation</li> <li>• Air pressure 700 hPa - 1060 hPa</li> </ul>
Transport and storage	<ul style="list-style-type: none"> <li>• -25°C to 70°C</li> <li>• Humidity up to 93%, without condensation</li> <li>• Air pressure 700 hPa - 1060 hPa</li> <li>• after external transport, you should wait 1-2 minutes until the device has reached room temperature</li> </ul>
Read the instruction manual.	
Date of manufacture	Month/year serial number

## 7. How to use the device

1. Switch on the device Press the I / GO button briefly.
2. Set the therapy time per STEP with the TIME button. Start with one minute per STEP, i.e. total therapy duration - 3 minutes. Increase by 1 minute per STEP for every subsequent therapy session. Do step 1 and step 2 (see below and following page) once a day.
3. Apply the device to the pain point.
4. Start therapy by pressing the I / GO button.
5. A beep sounds after each STEP.
6. The end of the therapy is signalled by 2 sounds.
7. Switch off the device: press and hold the I / GO button or switch off automatically after 30 seconds.

---

**Error message:**

***Bad or no skin contact is signalled with continuous beeps.***

---

### 1<sup>st</sup> Step:

Local therapy at the point of pain

## 2<sup>nd</sup> Step: Therapy in the abdominal area

Ensure good physical contact. If necessary, moisten with water. Drink plenty of clean water before and after the therapy.

## 8. Indications

- Pain therapy

## 9. Contraindications

- Pregnancy: Do not use the device.
- Patients with pacemakers: Do not use the device.
- Open wounds: Do not use the device.
- Fresh scars: Do not use the device.

# 10. Safety instructions for use

The following symbols may be given on the device and accessories:

## Symbols - Marking and their meaning

---

Symbol	Explanation
	Non-ionising electromagnetic radiation  <b>tested according to: IEC 60601-1-2:2014/AMD1:2020</b>
 <b>Medical device</b>	Indicates that the article is a medical device.
 <b>Serial number</b>	Indicates the manufacturer's serial number so that a specific medical device can be identified.
 <b>Keep dry</b>	Refers to a medical device that must be protected from moisture.
 <b>Manufacturer</b>	Indicates the manufacturer of the medical device.

- Please note that portable and mobile high frequency communications systems may cause interference to this device, even if they comply with the CISPR emission requirements. Portable high frequency communications equipment must not be used closer than 30 cm from the device.
- Using the POWERTUBE in the immediate vicinity of a shortwave or microwave device can lead to instability of the device output.
- The POWERTUBE may only be used for the electrical stimulation of nerves and muscles.
- The device is not designed for oxygen-rich environments.
- The POWERTUBE must not come into contact with flammable substances.
- Watch out for strangulation risks from the cable.
- Do not use on patients with pacemakers.
- Not suitable for patients under 2 years of age.
- Keep the therapy time short during the first use. (2-2-2 minutes)
- After one week, the therapy time can be increased slightly.
- Simultaneous connection with other devices is not permitted.
- POWERTUBE may only be operated with a 9 volt alkaline battery.
- POWERTUBE may only be operated with the original accessories included in the delivery.
- Dispose of POWERTUBE as hazardous waste.
- For technical questions, handling or maintenance of the POWERTUBE, please contact the manufacturer.

---

***Swisspowertube AG can be responsible for the reliability and safety of POWERTUBE only if:***

***“POWERTUBE is operated in accordance with the instruction manual,  
“extensions and modifications are solely carried out by persons authorised by Swisspowertube AG***

## **11. Care, maintenance and storage**

Effective protection of the POWERTUBE against damage includes:

- Proper operation and maintenance.
- Protecting the device from moisture, high temperatures, dust and direct sunlight .
- Keeping away from children, animals, and pests.

After every use, clean the contact surfaces of the POWERTUBE device and the hand electrodes with a commercially available cleaning agent (e.g. 1-propanil, ethanol, glyoxal). The printed type plate should not be cleaned.

- If the device is not used for a long period of time, remove the battery from the device.
- Dispose of batteries as hazardous waste.
- The expected service life of new alkaline batteries is approx. 2 months (with daily use for 20 minutes).

## 12. Malfunctions, function test

- When the battery runs out, this is indicated by a battery indicator LED.
- After switching on, POWERTUBE automatically does a battery test.
- If the hand electrode or the POWERTUBE body is not correctly positioned during therapy, a warning signal will sound.

### Do the function test yourself

1. Use a new battery for the self-test.
2. Connect the hand electrode to the cable.
3. Switch on the POWERTUBE.
4. Set the therapy time to 1 minute (default setting).
5. Start the therapy and make contact with the hand electrode on the POWERTUBE body (not on the therapy head) .
6. If you now hear a continuous beep, the POWERTUBE is functioning correctly.
7. Switch off POWERTUBE and restart therapy.

## 13. Accessories

- Hand electrode Gold (Art. No. 1018), silver (Art. No. 1019)
- Cable of the Powertube (Art. No. 1016)
- Instruction manual
- 9 volt alkaline battery 6LR61
- Softbag Gold (Art. No. 1010), silver (Art. No. 1011)

## 14. Guarantee

- The device and accessories must be inspected for defects and transport damage as soon as they are received. Claims for damages in this regard can only be asserted if the seller or forwarding agent is notified immediately. A damage report must be prepared.

- The warranty period for POWERTUBE is 2 years from the date of delivery.
- The scope of the guarantee extends to the free repair of the device, provided that the damage was caused by a material or manufacturing defect.
- Further claims, in particular claims for damages, are excluded.
- Repairs may only be carried out by Swisspowertube AG, its representatives or authorised specialist dealers.
- The warranty claim is void if improper modifications or repairs have been carried out.
- The warranty is neither extended nor renewed by services under warranty.
- Swisspowertube AG reserves all rights to this instruction manual, in particular the right of reproduction, distribution and translation. No part of this instruction manual may be reproduced or distributed in any form without the prior written consent of Swisspowertube AG.
- The information contained in this instruction manual may be changed or extended by Swisspowertube AG.
- Updates to the instruction manual are available for download on the website [www.swisspowertube.ch](http://www.swisspowertube.ch).

## 15. Safety check

In order to prevent a reduction in the safety of the device due to ageing, wear and tear, etc., section 6 of the Ordinance on the Installation, Operation and Use of Active Medical Devices prescribes regular safety checks.

The user must carry out the safety checks specified for this device in due time and to the prescribed extent. The safety checks may only be carried out by the manufacturer or by persons expressly authorised by the manufacturer.

## The following technical safety checks have been defined for POWERTUBE:

- Reading the labelling and the instruction manual
- Visual inspection of the device and accessories for damage
- Function check of the buttons and LED indicators

## 16. EC Declaration of Conformity CE 0482

Name and address of the manufacturer	Swisspowertube AG Breite 9 CH-3636 Forst (Switzerland) SRN no: CHRN-MF-20000238
As the manufacturer, we declare under our sole responsibility that	
the medical devices	POWERTUBE (gold) type FPT.B1 and POWERTUBE (silver) type FPT-D.B2
of the class	Ila Rule 9 according to Annex IX of Directive 93/42/EEC
comply with the relevant provisions of the Medical Devices Directive 93/42/EEC and its transposition into national law. The declaration of conformity refers to the medical devices supplied.	
Conformity assessment procedure	Directive 93/42/EEC Annex II, without section 4
Notified body	MEDCERT GmbH Pilatuspool 2 D-20355 Hamburg  <b>CE</b> 0482
CH-Forst, 24 <sup>th</sup> August 2021	Martin Frischknecht, CEO  
Name and address of the EU authorised representative	Wellbalance GmbH Klosterstrasse 10 D-77716 Haslach Germany SRN no: DE-AR-000010274

Technical description EN/IEC 60601-1-2


Table 1

Guideline and manufacturer's declaration - Electromagnetic emissions		
<p>The device is intended for use in the ELECTRONIC ENVIRONMENT specified below. The user should ensure that it is used in such an environment. The medical device runs under IEC 60601-1-2:2014/AMD:2020</p>		
Emission measurements	Compliance	Electromagnetic environment - Guidelines
High frequency emission according to CISPR 11	Group 1	The device uses high frequency energy exclusively for its internal FUNCTION. Therefore, its high frequency emission is very low and it is unlikely to interfere with neighbouring electronic devices.
High frequency emission according to CISPR 11	Class B	The device is intended for use in all establishments, including in residential areas and those directly connected to a PUBLIC POWER SUPPLY NETWORK that supplies buildings used for residential purposes.
Harmonics according to IEC 61000-3-2	n/a	Battery operated
Voltage fluctuations/flicker according to IEC 61000-3-3	n/a	

**Guideline and manufacturer's declaration -  
Electromagnetic immunity**

The device is intended for use in the ELECTRONIC ENVIRONMENT specified below. The user should ensure that it is used in such an environment.

Electrostatic discharge (ESD) IEC 61000-4-2	+/-8 kV contact, +/-15 kV air	
Emitted high frequency IEC 61000-4-3	80 MHz to 6 GHz	10 V/m
	380 MHz to 390 MHz	27 V/m
	430 MHz to 470 MHz	28 V/m
	704 MHz to 787 MHz	99 V/m
	800 MHz to 960 MHz	28 V/m
	1700 MHz to 1990 MHz	28 V/m
	2400 MHz to 2570 MHz	28 V/m
	5100 MHz to 5800 MHz	9 V/m
Conducted IEC 61000-4-6	0.15 MHz to 80 MHz	3 Vrms
ISM/Amateur bands	150 kHz to 80 MHz	6 Vrms
	Only for the electrode cable	
Mains frequency (50/60 Hz) Magnetic field IEC 61000-4-8	30A/m	
Magnetic proximity fields IEC 61000-4-39	<b>Test level [A/m]</b>	
<b>Frequencies</b>		
30 kHz CW	8	
134.2 kHz Pulse modulation 2.1 kHz	65	
13.56 MHz Pulse modulation 50 kHz	7.5	

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©2024 Swisspowertube AG. All rights reserved. All trade names are registered trademarks of

Swisspowertube AG  
Breite 9  
3636 Forst  
Switzerland

Edition 2024-01-26  
Printed in Switzerland

# Power+tube<sup>®</sup>